



**Lift Chair**

Useful for person after surgery, or with handicap or elderly who have difficulty standing.


STAND makes motor extend so chair mechanism moves from recline position through to standing.

SIT retracts motor and moves from stand, to sit, to recline position

STOP allows person to control exact position

**CABLE**

Important to use only this type of cable. It is the most common for USA chairs.



Back side as complete in the carton

Hook accessory  
Includes 1 plastic part and 1 screw

Clip accessory  
Includes plastic part top part with a spring, pivot pin, lower rotating plastic plate, and 1 screw



Finished Controller



Voice Controller in the Package